

# Cytoplast® PTFE Suture

The *soft* monofilament suture

- Size: USP 2-0
- 19 mm 3/8 Circle Reverse Cutting Needle  
(12 sutures per box)



**CS0418**

USP 2-0 with a 19 mm reverse cutting needle is popular for procedures where a longer needle is desired.

**Recommended**

- Size: USP 3-0
- 19 mm 3/8 Circle Reverse Cutting Needle  
(12 sutures per box)



**CS051819**

USP 3-0 with a 19 mm reverse cutting needle combines the most popular needle in dentistry with our most popular suture size.

- Size: USP 3-0
- 16 mm 3/8 Circle Reverse Cutting Needle  
(12 sutures per box)



**CS0518**

USP 3-0 with a 16 mm reverse cutting needle is by far the most popular size, especially for dental implant and bone grafting procedures.

- Size: USP 4-0
- 16 mm 3/8 Circle Reverse Cutting Needle  
(12 sutures per box)



**CS0618RC**

USP 4-0 with a 16 mm reverse cutting needle is the most popular of the smaller 4-0 sutures, especially for dental implant and bone grafting procedures.

- Size: USP 4-0
- 13 mm 1/2 Circle Taper Point Needle  
(12 sutures per box)



**CS0618PERIO**

USP 4-0 with a 13 mm taper point needle is often used on soft tissue grafts and suturing delicate tissues that require an atraumatic needle.

- Size: USP 4-0
- 13 mm 3/8 Circle Premium Reverse Cutting Needle  
(12 sutures per box)



**CS0618PREM**

USP 4-0 with a 13 mm premium reverse cutting needle offers over 50% longer cutting edges, giving needles a more slender, defined shape that greatly decreases the force applied to the needle during surgery.

# Cytoplast® PTFE Suture

The *soft* monofilament suture

## Features & Benefits

100% Medical Grade PTFE	Biologically inert
Monofilament	Doesn't wick bacteria
Soft (not stiff like most monofilaments)	Comfortable for patients
Little to no package memory	Excellent handling, knots securely
Non-resorbable	Keeps the surgical site reliably closed

## Praise for our Product

*"For years I have sworn by GORE-TEX® suture when I wanted insurance that a flap would not open during the healing period. Cytoplast® PTFE suture provides the same insurance at a better cost."*

*Cornell McCullom III, DDS, MD; Oral and Maxillofacial Surgeon*

*"Closing the surgery sites with Cytoplast® PTFE is easy and fast, and the patients hardly know there are sutures in the mouth."*

*Thomas Reinhart, DDS; Periodontist*

*"I love this stuff. It's easy to use, lasts as long as I need in the mouth, and doesn't wick infection."*

*Kurt Schneider, DDS; General Practitioner*

*"Keeps the surgical site reliably closed."*  
*Allen Schultz, DDS; Periodontist*

*GORE-TEX® is a registered trademark of W. L. Gore and Associates, Inc.*



Little to no package memory.



Excellent handling, knots securely.