



# OrACELL<sup>®</sup>

## Acellular Dermis

### Clinical Overview

An acellular dermal matrix used for soft tissue reconstruction in periodontal defects and ridge preservation

### Applications

- Guided Bone Regeneration: Extraction Socket, Ridge Augmentation, Sinus Augmentation
- Soft Tissue: Root Coverage, Insufficient Attached Gingiva

### Why Use

- **Safe** – Sterile, with a Sterility Assurance Level (SAL) of  $10^{-6}$
- **Decellularized** – Uses Matracell<sup>®</sup> technology to remove DNA and cellular remnants to decrease likelihood of an immune response
- **Biocompatible** –  $\geq 97\%$  of DNA removed and retains natural growth factors, collagen and elastin
- **No re-hydration needed** – Supplied pre-hydrated, ready-to-use
- **Convenient storage** – Room temperature (15°C to 30°C)





## Oracell Acellular Dermis

Thickness (mm)	Size (mm)	Shelf Life	Order Code
0.76 - 1.25	1.0 x 1.0	18 Months	OCELL 150
	1.0 x 4.0	18 Months	OCELL 151
	1.5 x 2.0	18 Months	OCELL 100
	2.0 x 4.0	3 Years	OCELL 101
1.26 - 1.75	1.0 x 1.0	18 Months	OCELL 250
	1.0 x 4.0	18 Months	OCELL 251
	1.5 x 2.0	18 Months	OCELL 200
	2.0 x 4.0	3 Years	OCELL 201

Instructions for Use available at [www.lifenethealth.org/IFU](http://www.lifenethealth.org/IFU)

68-60-084-02.03

LifeNet Health, the LifeNet Health logo, Oracell, and Matracell are registered trademarks of LifeNet Health.  
©2019 LifeNet Health, Virginia Beach, VA. All rights reserved.



**North America**  
1.888.847.7831  
[orders@lifenethealth.org](mailto:orders@lifenethealth.org)

**Europe**  
+ 43 1 375002710  
[eu\\_orders@lifenethealth.eu](mailto:eu_orders@lifenethealth.eu)

**Latin America** □ **Asia** □ **Middle East**  
1.757.464.4761 ext. 2000  
[internat.orders@lifenethealth.org](mailto:internat.orders@lifenethealth.org)

**LifeNetHealth.org**  
**LifeNetHealth.eu**